

DNAJC16 Antibody



PACO08888

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DNAJC16

Uniprot

Q9Y2G8

Synonyms:

DnaJ (Hsp40) homolog, subfamily C, member 16; DNAJC16; DKFZp686G1298; DKFZp686N0387; DKFZp781I1547; KIAA0962 ;

Immunogen:

Human DNAJC16.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A