

DNAJC27 Antibody



PACO08891

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DNAJC27

Uniprot

Q9NZQ0

Synonyms:

DnaJ (Hsp40) homolog, subfamily C, member 27; DNAJC27; DKFZp434N211; RBJ; RabJS ;

Immunogen:

Human DNAJC27.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A