

# DNAJC4 Antibody



PACO08892

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DNAJC4

**Uniprot**

Q9NNZ3

**Synonyms:**

DnaJ (Hsp40) homolog, subfamily C, member 4; DNAJC4; DANJC4; HSPF2; MCG18; MGC19482; MGC57189; MGC71863 ;

**Immunogen:**

Human DNAJC4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A