

# DNAJC5B Antibody



PACO08893

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DNAJC5B

**Uniprot**

Q9UF47

**Synonyms:**

DnaJ (Hsp40) homolog, subfamily C, member 5 beta; DNAJC5B; CSP-beta; MGC26226 ;

**Immunogen:**

Human DNAJC5B.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A