

DOCK9 Antibody



PACO08911

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

DOCK9

Uniprot

Q9BZ29

Synonyms:

dedicator of cytokinesis 9; DOCK9; DKFZp686C11110; DKFZp686D2047; DKFZp686N04132; FLJ11949; FLJ16744; FLJ44528; FLJ45282; FLJ45601; KIAA1058; KIAA1085; ZIZ1; ZIZIMIN1 ;

Immunogen:

Human DOCK9.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A