

DPYSL5 Antibody



PACO08932

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DPYSL5

Uniprot

Q9BPU6

Synonyms:

dihydropyrimidinase-like 5; DPYSL5; CRAM; CRMP-5; CRMP5; FLJ45383; Ulip6 ;

Immunogen:

Human DPYSL5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A