

DRD2 Antibody



PACO08936

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DRD2

Uniprot

P14416

Synonyms:

dopamine receptor D2; DRD2; D2DR; D2R ;

Immunogen:

Human DRD2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A