

DTX3L Antibody



PACO08947

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DTX3L

Uniprot

Q8TDB6

Synonyms:

deltex 3-like (Drosophila); DTX3L; BBAP ;

Immunogen:

Human DTX3L.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A