

DUS1L Antibody



PACO08949

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DUS1L

Uniprot

Q6P1R4

Synonyms:

dihydrouridine synthase 1-like (S. cerevisiae); DUS1L; DUS1; PP3111 ;

Immunogen:

Human DUS1L

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A