

DUS3L Antibody



PACO08950

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DUS3L

Uniprot

Q96G46

Synonyms:

dihydrouridine synthase 3-like (S. cerevisiae); DUS3L; DUS3; FLJ13896 ;

Immunogen:

Human DUS3L

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A