

DUSP10 Antibody



PACO08951

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DUSP10

Uniprot

Q9Y6W6

Synonyms:

dual specificity phosphatase 10; DUSP10; MKP-5; MKP5 ;

Immunogen:

Human DUSP10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A