

# DUSP15 Antibody



PACO08954

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

DUSP15

**Uniprot**

Q9H1R2

**Synonyms:**

dual specificity phosphatase 15; DUSP15; FLJ40111; VHY ;

**Immunogen:**

Human DUSP15.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A