

DUSP22 Antibody



PACO08956

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DUSP22

Uniprot

Q9NRW4

Synonyms:

dual specificity phosphatase 22; DUSP22; FLJ35864; JKAP; JSP1; MKPX; VHX ;

Immunogen:

Human DUSP22.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A