

DUSP4 Antibody



PACO08957

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DUSP4

Uniprot

Q13115

Synonyms:

dual specificity phosphatase 4; DUSP4; HVH2; MKP-2; MKP2; TYP ;

Immunogen:

Human DUSP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A