

DZIP1 Antibody



PACO08972

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DZIP1

Uniprot

Q86YF9

Synonyms:

DAZ interacting protein 1; DZIP1; DZIP; DZIPt1; KIAA0996; RP11-23E3.3 ;

Immunogen:

Human DZIP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A