

E2F4 Antibody



PACO08974

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

E2F4

Uniprot

Q16254

Synonyms:

E2F transcription factor 4, p107/p130-binding; E2F4; E2F-4 ;

Immunogen:

Human E2F4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A