

ECEL1 Antibody



PACO08984

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ECEL1

Uniprot

O95672

Synonyms:

endothelin converting enzyme-like 1; ECEL1; DINE; ECEX; XCE ;

Immunogen:

Human ECEL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A