

# EEF1A2 Antibody



PACO08998

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

EEF1A2

**Uniprot**

Q05639

**Synonyms:**

eukaryotic translation elongation factor 1 alpha 2; EEF1A2; EEF1AL; EF-1-alpha-2; EF1A; FLJ41696; HS1; STN; STNL ;

**Immunogen:**

Human EEF1A2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A