

PACO09011

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

EFS

Uniprot

O43281

Synonyms:

embryonal Fyn-associated substrate; EFS; CAS3; CASS3; EFS1; EFS2; HEFS; SIN ;

Immunogen:

Human EFS.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A