

EHBP1 Antibody



PACO09020

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

EHBP1

Uniprot

Q8NDI1

Synonyms:

EH domain binding protein 1; EHBP1; HPC12; KIAA0903; NACSIN ;

Immunogen:

Human EHBP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A