

EIF1AD Antibody



PACO09027

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

EIF1AD

Uniprot

Q8N9N8

Synonyms:

eukaryotic translation initiation factor 1A domain containing; EIF1AD; MGC11102 ; Haponin

Immunogen:

Human EIF1AD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A