

EIF2B4 Antibody



PACO09038

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

EIF2B4

Uniprot

Q9UI10

Synonyms:

eukaryotic translation initiation factor 2B, subunit 4 delta, 67kDa; EIF2B4; DKFZp586J0119; EIF-2B; EIF2B; EIF2Bdelta ;

Immunogen:

Human EIF2B4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A