

EIF3D Antibody



PACO09044

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

EIF3D

Uniprot

O15371

Synonyms:

eukaryotic translation initiation factor 3, subunit D; EIF3D; EIF3S7; MGC126526; MGC17258; eIF3-p66; eIF3-zeta ;

Immunogen:

Human EIF3D.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A