

EIF3E Antibody



PACO09045

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

EIF3E

Uniprot

P60228

Synonyms:

eukaryotic translation initiation factor 3, subunit E; EIF3E; EIF3-P48; EIF3S6; INT6; eIF3-p46 ;

Immunogen:

Human EIF3E.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A