

EIF4A3 Antibody



PACO09052

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

EIF4A3

Uniprot

P38919

Synonyms:

eukaryotic translation initiation factor 4A, isoform 3; EIF4A3; DDX48; KIAA0111; MGC10862; NMP265; NUK-34; eIF4AIII; hNMP265 ;

Immunogen:

Human EIF4A3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A