

EIF4EBP2 Antibody



PACO09056

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

EIF4EBP2

Uniprot

Q13542

Synonyms:

eukaryotic translation initiation factor 4E binding protein 2; EIF4EBP2;4EBP2; PHASII ;

Immunogen:

Human EIF4EBP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A