

EIF5A2 Antibody



PACO09062

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

EIF5A2

Uniprot

Q9GZV4

Synonyms:

eukaryotic translation initiation factor 5A2; EIF5A2; EIF-5A2; eIF5AII ;

Immunogen:

Human EIF5A2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A