

ELMOD3 Antibody



PACO09077

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ELMOD3

Uniprot

Q96FG2

Synonyms:

ELMO/CED-12 domain containing 3; ELMOD3; FLJ21977; FLJ35601; MGC111036; RBED1; RBM29 ;

Immunogen:

Human ELMOD3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A