

ELOVL2 Antibody



PACO09079

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ELOVL2

Uniprot

Q9NXB9

Synonyms:

elongation of very long chain fatty acid, (FEN1/Elo2, SUR4/Elo3, yeast)-like 2; ELOVL2; FLJ20334; SSC2 ;

Immunogen:

Human ELOVL2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A