

# ELP4 Antibody



PACO09082

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ELP4

**Uniprot**

Q96EB1

**Synonyms:**

elongation protein 4 homolog (S. cerevisiae); ELP4; C11orf19; FLJ20498; PAX6NEB; PAXNEB; dJ68P15A.1 ;

**Immunogen:**

Human ELP4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A