

# EML1 Antibody



PACO09087

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

EML1

**Uniprot**

O00423

**Synonyms:**

echinoderm microtubule associated protein like 1; EML1; ELP79; EMAP; EMAPL; FLJ45033; HuEMAP ;

**Immunogen:**

Human EML1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A