

# ENOPH1 Antibody



PACO09093

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ENOPH1

**Uniprot**

Q9UHY7

**Synonyms:**

enolase-phosphatase 1; ENOPH1; DKFZp586M0524; E1; FLJ12594; MASA; MST145 ;

**Immunogen:**

Human ENOPH1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A