

ENTPD8 Antibody



PACO09105

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ENTPD8

Uniprot

Q5MY95

Synonyms:

ectonucleoside triphosphate diphosphohydrolase 8; ENTPD8; GLSR2492; NTPDase-8; UNQ2492 ; NTPDase8

Immunogen:

Human ENTPD8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A