

EPB41L2 Antibody



PACO09107

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

EPB41L2

Uniprot

O43491

Synonyms:

erythrocyte membrane protein band 4.1-like 2; EPB41L2;4.1-G; DKFZp781D1972; DKFZp781H1755 ;

Immunogen:

Human EPB41L2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A