

# EPSTI1 Antibody



PACO09118

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

EPSTI1

**Uniprot**

Q96J88

**Synonyms:**

epithelial stromal interaction 1 (breast); EPSTI1; BRESI1; MGC29634 ;

**Immunogen:**

Human EPSTI1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A