

ERC2 Antibody



PACO09126

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ERC2

Uniprot

O15083

Synonyms:

ELKS/RAB6-interacting/CAST family member 2; ERC2; CAST; CAST1; ELKSL; KIAA0378; MGC133063; MGC133064; SPBC110; Spc110 ;

Immunogen:

Human ERC2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A