

ERGIC1 Antibody



PACO09134

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ERGIC1

Uniprot

Q969X5

Synonyms:

endoplasmic reticulum-golgi intermediate compartment (ERGIC) 1; ERGIC1; ERGIC-32; ERGIC32; FLJ39864; KIAA1181; MGC14345 ;

Immunogen:

Human ERGIC1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A