

ERGIC2 Antibody



PACO09135

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ERGIC2

Uniprot

Q96RQ1

Synonyms:

ERGIC and golgi 2; ERGIC2; CDA14; Erv41; MGC111152; PTX1; cd002 ;

Immunogen:

Human ERGIC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A