

ERGIC3 Antibody



PACO09136

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ERGIC3

Uniprot

Q9Y282

Synonyms:

ERGIC and golgi 3; ERGIC3; C20orf47; CGI-54; Erv46; NY-BR-84; PRO0989; SDBCAG84; dJ477O4.2;

Immunogen:

Human ERGIC3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A