

# ESX1 Antibody



PACO09149

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ESX1

**Uniprot**

Q8N693

**Synonyms:**

ESX homeobox 1; ESX1; ESX1L; ESXR1 ;

**Immunogen:**

Human ESX1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A