

# EXOC5 Antibody



PACO09164

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

EXOC5

**Uniprot**

O00471

**Synonyms:**

exocyst complex component 5; EXOC5; DKFZp666H126; HSEC10; PRO1912; SEC10; SEC10L1; SEC10P ;

**Immunogen:**

Human EXOC5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A