

# EXOC6 Antibody



PACO09165

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

EXOC6

**Uniprot**

Q8TAG9

**Synonyms:**

exocyst complex component 6; EXOC6; DKFZp761I2124; EXOC6A; FLJ1125; FLJ11251; MGC33397; SEC15; SEC15L; SEC15L1; SEC15L3; Sec15p ;

**Immunogen:**

Human EXOC6.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A