

FADS1 Antibody



PACO09191

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FADS1

Uniprot

O60427

Synonyms:

fatty acid, desaturase 1; FADS1; D5D; FADS6; FADSD5; FLJ38956; FLJ90273; LLCDL1; TU12 ;

Immunogen:

Human FADS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A