

FAM101A Antibody



PACO09197

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RFLNA

Uniprot

Q6ZTI6

Synonyms:

family with sequence similarity 101, member A; FAM101A; FLJ44614 ;

Immunogen:

Human FAM101A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A