

FAM118B Antibody



PACO09201

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FAM118B

Uniprot

Q9BPY3

Synonyms:

family with sequence similarity 118, member B; FAM118B; FLJ21103 ;

Immunogen:

Human FAM118B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A