

FAM124B Antibody



PACO09204

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FAM124B

Uniprot

Q9H5Z6

Synonyms:

family with sequence similarity 124B; FAM124B; FLJ22746 ;

Immunogen:

Human FAM124B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A