

FAM127A Antibody



PACO09205

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RTL8C

Uniprot

O15255

Synonyms:

family with sequence similarity 127, member A; FAM127A; CXX1; MAR8C; MART8C; MGC117411; Mar8; Mart8 ;

Immunogen:

Human FAM127A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A