

FAM156A Antibody



PACO09208

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FAM156A

Uniprot

Q8NDB6

Synonyms:

family with sequence similarity 156, member A; FAM156A; DKFZp686B22211; PRO0659; TMEM29 ;

Immunogen:

Human FAM156A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A