

FAM160B2 Antibody



PACO09209

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FAM160B2

Uniprot

Q86V87

Synonyms:

family with sequence similarity 160, member B2; FAM160B2; FLJ11125; FLJ21801; MGC138352; RAI16 ;

Immunogen:

Human FAM160B2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A