

FAM185A Antibody



PACO09211

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FAM185A

Uniprot

Q8N0U4

Synonyms:

family with sequence similarity 185, member A; FAM185A; MGC35361 ;

Immunogen:

Human FAM185A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A